

Item No.	Description
1	2
1	LP (SPINAL) NEEDLE
	Disposable, sterile pack
	Luer lock hub transperant
	hub lumber puncture
	needle. Length Range,
	8.5cm to 10cm, small
	bevel indicator for bevel
	direction.
	Guage
	23
2	Sterile disposable
	Hypodermic plastic
	syringe 2cc
	With Leur Lock
3	Sterile disposable
	Hypodermic plastic
	syringe 2cc
	Without Leur Lock

4	Sterile disposable
	Hypodermic plastic
	syringe 2cc
	With Needle.
5	Sterile Non-reusable Auto
	Disable Mechanism
	incorporated Syringes
	2 ml.
6	Sterile disposable
	Hypodermic plastic
	syringe 5cc, cn.
	With Leur Lock
7	Sterile disposable
	Hypodermic plastic
	syringe 5cc, cn.
	With Needle
8	Sterile Non-reusable Auto
	Disable Mechanism
	incorporated Syringes
	5 ml.

9	Sterile disposable
	Hypodermic plastic
	syringe 10cc, sn.
	With Leur Lock
10	Sterile disposable
	Hypodermic plastic
	syringe 10cc, sn.
	With Needle.
11	Sterile Non-reusable Auto
	Disable Mechanism
	incorporated Syringes
	10 ml.
12	Sterile disposable
	Hypodermic plastic
	syringe 20 cc, sn.
	without Leur Lock
13	Sterile disposable
	Hypodermic plastic
	syringe 20 cc, sn.
	with Needle, size 20 G,

	1 1/2" long.
14	Sterile disposable
	Hypodermic plastic
	syringe 50 cc, sn.
	With Leur Lock
15	Tuberculin Syringes,
	Disposable sterile, sterile
	pack, Leur Lock, hub 1ml
16	Disposable needle -
	Sterile size - 18, length
	1 1/2"
17	Disposable needle -
	Sterile size - 20, length
	1 1/2"
18	Disposable needle sterile
	Length 1 1/2"
	Gauge - 21
19	Epidural Catheters with
	Filter. 0.2mm, flat Sterile
	Pack with marking clear

	cathetar with lateral eyes
	at the tip, Luer Lock twist
	connector Cathetar
	guide/introducer
	Guage 16 - 25%
20	IV Cannula Disposable
	Sterile pack, Fluro
	Ethylene Polypropylene/or
	PEU vialon Radiopaque
	Cathetar with taper tip
	Distance between bevel
	of stilette and cannula
	should be minimal
	preferable to have flanges
	with no Injection ports.
	Guage - 14
	Length - 45mm to 75mm
21	IV Cannula Disposable
	Sterile pack, Fluro
	Ethylene Polypropolene/

	Radiopaque Cathetar
	with taper tip. Distance
	between bevel of stilette
	and cannula should be
	minimal preferable to
	have flanges with no
	Injection ports.
	Guage - 20
	Length - 25mm to 50mm
22	IV Cannula Disposable
	Sterile Fluro Ethylene
	Polypropolene/
	Radiopaque Cathetar with
	taper tip. Distance
	between bevel of stilette
	and cannula should be
	minimal preferable to
	have flanges with no
	Injection ports.
	Guage - 24

	Length - 18mm - 30mm
23	Complete Epidural Set
	(Disposable). Epidural
	needle, catheter filter &
	Syringe, Tuohy needle,
	Catheter with guide &
	lateral eyes marking
	gradation at 50-150mm
	at 10mm & at 200mm.
	Tuohy needle marking at
	3cms from tip & later on
	every 1 cm. Filter 0.2 mm,
	flat Loss of resistance
	syringe 10cc.
	Guage Length
	16 75mm - 110mm
24	Combined Spinal Epidural
	Needles Spinal epidural
	set pencil point/Lancet
	Disposable Spinal needle,

	sterile pack, 26G/16G,
	Luer lock & transparent
	hub, Tuohy needle
	length = 7.5 to 11cms.
	Spinal needle protrusion
	by 11 - 14mm with
	epidural Catheter, filter,
	Loss of resistance
	syringe.
25	Triple Lumen Central
	Venous Catheterisation
	Set for Adults
	Containing
	Catheter Size 7 Fr.
	3 Lumens : 16G, 18G,
	Length - 14 to 18 Cm.
	Latex free, Radio Opaque
	catheter, clamp / Extension
	line clamp 18G Needle,
	J Tip guide wire,

	Introducer - 1
	Tissue Dilator - 1
	Syringe : Leur Lock /
	Leur Slip Suture Wing
	Sterile, Non-pyrogenic,
	single use, peel to open
	pack, transperant cover
	on one side
26	Triple Lumen Central
	Venous Catheterisation
	Set for Children.
	Containing
	Catheter Size 5.5 Fr.
	3 Lumens : 20G, 22G,
	Length - 8 to 10 Cm.
	Latex free, Radio Opaque
	catheter / Extension
	line clamp, introducer
	Needle 21G,
	J Tip guide wire,

	Tissue Dilator - 1
	Syringe : Leur Lock /
	Leur Slip Suture Wing.
	Sterile, Non-pyrogenic,
	single use, peel to open
	pack, transperant cover
	on one side
27	LP (SPINAL) NEEDLE
	Disposable, sterile pack
	Luer lock hub transperant
	hub lumber puncture
	needle. Length Range,
	8.5cm to 10cm, small
	bevel indicator for bevel
	direction.
	Guage
	25
28	LP (SPINAL) NEEDLE
	Disposable, sterile pack
	Luer lock hub transperant

	hub lumbar puncture
	needle. Length Range,
	8.5cm to 10cm, small
	bevel indicator for bevel
	direction.
	Guage
	26
29	Epidural Catheters with
	Filter. 0.2mm, flat Sterile
	Pack with marking clear
	cathetar with lateral eyes
	at the tip, Luer Lock twist
	connector Cathetar
	guide/introducer
	Guage 18 - 50%
30	IV Cannula Disposable
	Sterile pack, Fluro
	Ethylene Polypropylene/or
	PEU vialon Radiopaque
	Cathetar with taper tip

	Distance between bevel
	of stilette and cannula
	should be minimal
	preferable to have flanges
	with no Injection ports.
	Guage - 16
	Length - 45mm to 75mm
31	IV Cannula Disposable
	Sterile pack, Fluro
	Ethylene Polypropylene/or
	PEU vialon Radiopaque
	Cathetar with taper tip
	Distance between bevel
	of stilette and cannula
	should be minimal
	preferable to have flanges
	with no Injection ports.
	Guage - 18
	Length - 45mm to 75mm
32	IV Cannula Disposable

	Sterile pack, Fluro
	Ethylene Polypropolene/
	Radiopaque Cathetar
	with taper tip. Distance
	between bevel of stilette
	and cannula should be
	minimal preferable to
	have flanges with no
	Injection ports.
	Guage - 22
	Length - 25mm to 50mm
33	Peripherally placed
	Radiopaque central
	Catheter. Guage 12/14
	(5%), 14/16 (50%),
	16/18 (20%),18/20
	(25%) (Needles Size/
	Catheter Size)
	Length - 70cm - 70%
34	Peripherally placed

	Radiopaque central
	Catheter. Gauge 12/14
	(5%), 14/16 (50%),
	16/18 (20%), 18/20
	(25%) (Needles Size/
	Catheter Size)
	Length - 45cm - 15%
35	Complete Epidural Set
	(Disposable). Epidural
	needle, catheter filter &
	Syringe, Tuohy needle,
	Catheter with guide &
	lateral eyes marking
	gradation at 50-150mm
	at 10mm & at 200mm.
	Tuohy needle marking at
	3cms from tip & later on
	every 1 cm. Filter 0.2 mm,
	flat Loss of resistance
	syringe 10cc.

	Guage	Length
	18	75mm - 110mm